

# EZEKLEEN 2.5 HD

4 - STAGE ULTRA CLEAN WATER SYSTEM

Model: 8000

## CHANGING YOUR CARTRIDGES - INSTRUCTIONS

### CHANGING CARTRIDGES BEFORE CHANGING 8002-FR, 8003 or 8004 CARTRIDGE

- A.1 Have a towel available to absorb any water released by system
- A.2 Have full tank of water before changing any cartridges
- A.3 Turn off water supply under sink by closing Ezee Kleen Shut Off Valve or cold water shut off valve
- A.4 Close Ezee Kleen Tank Valve
- A.5 Open Ezee Kleen Dispensing Faucet
- A.6 This will relieve line pressure in the system
- A.7 After completing any cartridge change wait for tank to fill and, with faucet closed, system will be fully pressurized, check for leaks

### CHANGE 8003 PRE FILTER CARTRIDGE ANNUALLY

- B.1 Remove Slip-in-Elbow from both ends of cartridge - Consult User's Guide to understand quick connect feature on all Ezee Kleen fittings
- B.2 Remove protective dust plugs from new cartridge
- B.3 Replace expired cartridge ensuring water is flowing in the correct direction
- B.4 Reinsert Slip-in-Elbows to both ends of cartridge
- B.5 Turn on water supply
- B.6 Open Ezee Kleen Tank Valve
- B.7 Run water from Ezee Kleen Dispensing Faucet until you have a steady stream of water - Continue until tank is empty
- B.8 Close Dispensing Faucet
- B.9 Allow 2 hours for tank to fill again - System is now ready to use

### CHANGE 8004 REVERSE OSMOSIS (RO) MEMBRANE BI-ANNUALLY (Every 2 years)

*(Note: Do not replace 8002-FR DI Cartridge at same time as replacing 8004 RO Membrane. Empty, and refill one full tank of water through new RO Membrane prior to changing 8002-FR DI cartridge)*

- C.1 Disconnect waterlines from all 3 fittings at the ends of RO Cartridge - Note which waterline each fitting is attached to
- C.2 Remove RO Cartridge from clips
- C.3 Unscrew cap of RO Cartridge and discard expired RO Membrane
- C.4 Clean inside of RO Cartridge with mild detergent - Rinse thoroughly - Insert new RO Membrane
- C.5 Screw on RO Cartridge cap - Do not over tighten - Fit RO Cartridge back into clips
- C.6 Reconnect all 3 waterlines to fittings on RO Cartridge
- C.7 Turn on water supply
- C.8 Open Ezee Kleen Tank Valve
- C.9 Run water from Ezee Kleen Dispensing Faucet until you have a steady stream of water - Continue until tank is empty
- C.10 Close Dispensing Faucet
- C.11 Allow 2 hours for tank to fill again - System is now ready to use

### CHANGE 8005 ULTRA-VIOLET (UV) LAMP ANNUALLY

*(Note: UV Lamp can be changed without turning off water)*

- D.1 Disconnect power supply to UV Cartridge
- D.2 Swivel/swing UV Cartridge outwards away from unit
- D.3 Grasp white electrical cable exiting end of UV Cartridge
- D.4 Gently pull electrical cable out of UV Cartridge - The entire lamp is attached to the electrical cable and will slide out of UV Cartridge
- D.5 Replace UV Lamp by sliding new UV Lamp into the UV Cartridge
- D.6 Re-connect power supply to UV Cartridge

### CHANGE 8002-FR DI CARTRIDGE IF MONITOR IS FLASHING WHILE WATER IS BEING DISPENSED FROM EZEK KLEEN DISPENSING FAUCET

**Monitor will also be reading >5.00 - Above this level water is unacceptable for autoclave use**

*(Note: Monitor can accurately gauge water quality to 19.99. Above 19.99 monitor will show "1 ." and is still unacceptable for autoclave use)*

- E.1 Remove Slip-in-Elbow from both ends of cartridge - Consult User's Guide to understand quick connect feature on all Ezee Kleen fittings
- E.2 Remove protective dust plugs from new cartridge
- E.3 Replace expired cartridge ensuring water is flowing in the correct direction
- E.4 Reinsert Slip-in-Elbows to both ends of cartridge
- E.5 Replace 8006 Drainline Flow Restrictor on black 1/4" drainline, ensuring arrow is pointing away from the unit
- E.6 Turn on water supply
- E.7 Open Ezee Kleen Tank Valve
- E.8 Run water from Ezee Kleen Dispensing Faucet until you have a steady stream of water - Continue until tank is empty
- E.9 Check that monitor is < 5.00. Monitor will stop flashing red
- E.10 Allow 2 hours for tank to fill again- System is now ready to use

Call Oasis Dental Group Inc. @ (800) 338-6693, M-F, 8:00 - 5:00 PST for technical assistance.